

# Technical Information Report



## AAMI TIR37: 2013

Sterilization of health care  
products—Radiation—  
Guidance on sterilization of  
biologics and tissue-based  
products



# **Sterilization of health care products— Radiation—Guidance on sterilization of biologics and tissue-based products**

Approved 30 December 2013 by  
**Association for the Advancement of Medical Instrumentation**

**Abstract:** This AAMI Technical Information Report (TIR) provides guidance for development, validation and routine control associated with the radiation sterilization processing of biologics and tissue-based products.

**Keywords:** sterilization, radiation, tissue

## AAMI Technical Information Report

A technical information report (TIR) is a publication of the Association for the Advancement of Medical Instrumentation (AAMI) Standards Board that addresses a particular aspect of medical technology.

Although the material presented in a TIR may need further evaluation by experts, releasing the information is valuable because the industry and the professions have an immediate need for it.

A TIR differs markedly from a standard or recommended practice, and readers should understand the differences between these documents.

Standards and recommended practices are subject to a formal process of committee approval, public review, and resolution of all comments. This process of consensus is supervised by the AAMI Standards Board and, in the case of American National Standards, by the American National Standards Institute.

A TIR is not subject to the same formal approval process as a standard. However, a TIR is approved for distribution by a technical committee and the AAMI Standards Board.

Another difference is that, although both standards and TIRs are periodically reviewed, a standard must be acted on—reaffirmed, revised, or withdrawn—and the action formally approved usually every five years but at least every 10 years. For a TIR, AAMI consults with a technical committee about five years after the publication date (and periodically thereafter) for guidance on whether the document is still useful—that is, to check that the information is relevant or of historical value. If the information is not useful, the TIR is removed from circulation.

A TIR may be developed because it is more responsive to underlying safety or performance issues than a standard or recommended practice, or because achieving consensus is extremely difficult or unlikely. Unlike a standard, a TIR permits the inclusion of differing viewpoints on technical issues.

**CAUTION NOTICE:** This AAMI TIR may be revised or withdrawn at any time. Because it addresses a rapidly evolving field or technology, readers are cautioned to ensure that they have also considered information that may be more recent than this document.

All standards, recommended practices, technical information reports, and other types of technical documents developed by AAMI are *voluntary*, and their application is solely within the discretion and professional judgment of the user of the document. Occasionally, voluntary technical documents are adopted by government regulatory agencies or procurement authorities, in which case the adopting agency is responsible for enforcement of its rules and regulations.

Comments on this technical information report are invited and should be sent to AAMI, Attn: Standards Department, 4301 N. Fairfax Drive, Suite 301, Arlington, VA 22203-1633.

### *Published by*

Association for the Advancement of Medical Instrumentation  
4301 N. Fairfax Drive, Suite 301  
Arlington, VA 22203-1633  
www.aami.org

© 2014 by the Association for the Advancement of Medical Instrumentation

All Rights Reserved

Publication, reproduction, photocopying, storage, or transmission, electronically or otherwise, of all or any part of this document without the prior written permission of the Association for the Advancement of Medical Instrumentation is strictly prohibited by law. It is illegal under federal law (17 U.S.C. § 101, *et seq.*) to make copies of all or any part of this document (whether internally or externally) without the prior written permission of the Association for the Advancement of Medical Instrumentation. Violators risk legal action, including civil and criminal penalties, and damages of \$100,000 per offense. For permission regarding the use of all or any part of this document, complete the reprint request form at [www.aami.org](http://www.aami.org) or contact AAMI at 4301 N. Fairfax Drive, Suite 301, Arlington, VA 22203-1633. Phone: +1-703-525-4890; Fax: +1-703-525-1067.

**ISBN 1–57020–522–1**

# Contents

Page

Glossary of equivalent standards.....	iv
Committee representation.....	v
Foreword.....	viii
Introduction.....	ix
1 Scope.....	1
1.1 Inclusions.....	1
1.2 Exclusions.....	1
2 Normative references.....	1
3 Terms and definitions.....	1
4 Definition and maintenance of product families for dose setting, dose substantiation and sterilization dose auditing.....	3
4.1 General.....	3
4.2 Defining product families.....	4
4.3 Designation of product to represent a product family for performance of a verification dose experiment or a sterilization dose audit.....	4
4.4 Maintaining product families.....	5
4.5 Effect of failure on establishment of sterilization dose or of a sterilization dose audit on a product family.....	5
5 Selection and testing of product for establishing and verifying the sterilization dose.....	5
5.1 General.....	5
5.2 Nature of product.....	6
5.3 Sample Item Portion (SIP).....	6
5.4 Sampling conditions.....	6
5.5 Microbiological testing.....	7
6 Methods of dose establishment.....	13
6.1 General.....	13
6.2 Selecting a method for the establishment of a sterilization dose.....	13
6.3 Establishing a maximum acceptable dose.....	13
Selection of a Dose Setting Method.....	14
Bibliography.....	16

## **Glossary of equivalent standards**

International Standards adopted in the United States may include normative references to other International Standards. AAMI maintains a current list of each International Standard that has been adopted by AAMI (and ANSI). Available on the AAMI website at the address below, this list gives the corresponding U.S. designation and level of equivalency to the International Standard.

[www.aami.org/standards/glossary.pdf](http://www.aami.org/standards/glossary.pdf)

## Committee representation

### Association for the Advancement of Medical Instrumentation Radiation Sterilization Working Group

This Technical Information Report (TIR) was developed by the AAMI Radiation Sterilization Working Group under the auspices of the AAMI Sterilization Standards Committee. Working Group approval of the TIR does not necessarily imply that all committee members voted for its approval.

At the time this document was published, the **AAMI Radiation Sterilization Working Group** had the following members:

**Chairs:** Trabue D. Bryans, BryKor LLC  
Pat Weikel, FDA/CDRH

**Members:** Edward Arscott, NAMSA  
Simon Bogdanský, PhD, Allo Source  
Curt Andrew Bogue, BAS, Cook Inc.  
Anne F. Booth, MS, Booth Scientific Inc.  
Harry F. Bushar, PhD  
David Cardin, Zimmer Inc.  
Denise Cleghorn, Boston Scientific Corporation  
Gary N. Cranston, Consulting & Technical Services/PCS  
Emily Craven, Nordion Inc.  
Greg Crego, Moog Medical Devices  
Elaine Daniell, CR Bard  
Douglas D. Davie, Sterilization Validation Services  
Darci Diage, Direct Flow Medical Inc.  
April J. Doering, St Jude Medical Inc.  
Niki Fidopiastis, Sterigenics International  
William F. FitzGerald, PE, FitzGerald & Associates Ltd  
Lisa Foster, Medpoint LLC  
Sarah Gagnon, Microtest Laboratories Inc.  
Shelley Green, WuXi AppTec Inc.  
Joyce M. Hansen, Johnson & Johnson  
Thomas L. Hansen, Terumo Americas Corporate  
Douglas F. Harbrecht, Sterility Assurance LLC  
Deborah A. Havlik, Hospira Worldwide Inc.  
Donna Horner, Abbott Laboratories  
Betty Howard, MBA MS BS, Steris Corporation  
Jim Kaiser, Bausch & Lomb Inc.  
David King, Tandem Diabetes Care Inc.  
Carolyn L. Kinsley, LexaMed Ltd  
Reynaldo Lopez, Cardinal Health (MP&S)  
Antonio Lopez-Feliciano, CareFusion  
Ronald G. Lulich, 3M Healthcare  
Jeff Martin, Alcon Laboratories Inc.  
Astrid Merrifield, Medtronic Inc. WHQ Campus  
Russell D. Mills, GE Healthcare  
Larry Nichols, Nutek Corporation  
Gerry A. O'Dell, MS, Gerry O'Dell Consulting  
Dave Parente, Ecolab  
Michelle Peterson, Stryker Instruments Division  
Rudy M. Pina, Dynatec Scientific Labs Inc.  
Jody Rupert, WL Gore & Associates Inc.  
Manuel Saavedra, Jr., Kimberly-Clark Corporation  
Jon Seulean, Terumo BCT  
Harry L. Shaffer, Sterilization Consulting Services

Michael G. Sprague, Ethide Laboratories Inc.  
Sopheak Srun, MPH SM (NRCM), Quality Tech Services Inc.  
Fenil Sutaria, Medline Industries Inc.  
Richard L. Weisman, Fresenius Medical Care Renal Therapies Group  
Patrick B. Weixel, FDA/CDRH  
Beverly G. Whitaker, CQA RAB MBA, Indigo Consulting Group Inc.  
John A. Williams, Baxter Healthcare Corporation  
Martell Kress Winters, BS SM, Nelson Laboratories Inc.  
Cheryl Work, Becton Dickinson & Company

**Alternates:** Marjean Boyter, Fresenius Medical Care Renal Therapies Group  
Carolyn Braithwaite-Nelson, Spectranetics Corporation  
Susan Butler, Boston Scientific Corporation  
Glenn W. Calvert, BS, Johnson & Johnson  
Claudia Camp, Stryker Instruments Division  
Rebecca Campbell, Sterigenics International  
Jessica Desmond, Microtest Laboratories Inc.  
John DiCaro, Covidien  
Dave Dion, Cardinal Health (MP&S)  
Brian R. Drumheller, CR Bard  
Henry Hart, Zimmer Inc.  
Brent Huberty, Medtronic Inc. WHQ Campus  
Nichole Jackson, Ecolab  
Chris Johnson, Steris Corporation  
Chris Kobus, GE Healthcare  
Ezra Koski, Terumo BCT  
Vu Le, Abbott Laboratories  
Jody L. Markling, St Jude Medical Inc.  
Consuelo McChesney, Alcon Laboratories Inc.  
Patrick J. McCormick, PhD, Bausch & Lomb Inc.  
Joseph Mello, Ethide Laboratories Inc.  
Kimberly Patton, Becton Dickinson & Company  
Michelle Pierce, NAMSA  
Nancy Rakiewicz, Moog Medical Devices  
Robert R. Reich, BS MS, LexaMed Ltd  
Reginald Roberts, Medline Industries Inc.  
Mark Seybold, Baxter Healthcare Corporation  
Wendy Wangsgard, PhD, Nelson Laboratories Inc.  
Keisha Weaver, Kimberly-Clark Corporation  
Casimir John Woss, PhD, FitzGerald & Associates Ltd

At the time this document was published, the **AAMI Sterilization Standards Committee** had the following members:

**Chairs:** Victoria M. Hitchins, PhD, FDA/CDRH  
Michael H. Scholla, PhD, Dupont Protection Technologies

**Members:** Trabue D. Bryans, BryKor LLC  
Susan Butler, Boston Scientific Corporation  
Nancy Chobin, RN CSPDM, St Barnabas Healthcare System  
Charles Cogdill, Covidien  
Ramona Conner, RN MSN CNOR, Association of Perioperative Registered Nurses  
Jacqueline Daley, Sinai Hospital of Baltimore  
Kimbrell Darnell, CR Bard  
Gordon M. Ely, WuXi AppTec Inc.  
Lisa Foster, Medpoint LLC  
Joel R. Gorski, PhD, NAMSA  
Joyce M. Hansen, Johnson & Johnson  
Douglas F. Harbrecht, Sterility Assurance LLC  
Deborah A. Havlik, Hospira Worldwide Inc.  
Susan G. Klacik, CCSMC FCS ACE, IAHCSMM  
Byron J. Lambert, PhD, Abbott Laboratories



Colleen Patricia Landers, RN, Timmins & District Hospital  
Reynaldo Lopez, Cardinal Health (MP&S)  
Jeff Martin, Alcon Laboratories Inc.  
Patrick J. McCormick, PhD, Bausch & Lomb Inc.  
Gerald E. McDonnell, PhD, Steris Corporation  
Janet M. Prust, 3M Healthcare  
Nancy Rakiewicz, Moog Medical Devices  
Mark Seybold, Baxter Healthcare Corporation  
Andrew Sharavara, PhD, Propper Manufacturing Co Inc.  
Mark N. Smith, Getinge USA  
James Sidney Wiggs, BSN CRCST, Legacy Health System  
Martell Kress Winters, BS SM, Nelson Laboratories Inc.  
William T. Young, Sterigenics International

**Alternates:** Lloyd Brown, Covidien  
Peter A. Burke, PhD, Steris Corporation  
Dave Dion, Cardinal Health (MP&S)  
Ken Eddington, NAMSA  
Thomas J. Frazar, Johnson & Johnson  
Martha M. Kadas, Sterigenics International  
Jim Kaiser, Bausch & Lomb Inc.  
Stacy Kromenhoek, Boston Scientific Corporation  
Natalie Lind, IAHCSSM  
Mary S. Mayo, CR Bard  
David Ford McGoldrick, BS, Abbott Laboratories  
Jerry R. Nelson, PhD, Nelson Laboratories Inc.  
Patrick Polito, Moog Medical Devices  
Karen Polkinghorne, Dupont Protection Technologies  
Mike Sadowski, Baxter Healthcare Corporation  
Craig A. Wallace, 3M Healthcare

---

NOTE—Participation by federal agency representatives in the development of this document does not constitute endorsement by the federal government or any of its agencies.

---

## Foreword

The dose setting methods described in the ANSI/AAMI/ISO 11137 series *Sterilization of health care products—Radiation*, Parts 1–3, were developed in the context of medical devices, and do not address the unique issues associated with biologics/tissues. These unique issues might potentially involve every aspect of a radiation sterilization validation, routine processing and maintenance of the sterilization process. This TIR addresses these unique issues and provides guidance on how to adapt current guidelines and standards for use with biologics/tissues.

As used within the context of this document, “should” indicates that among several possibilities, one is recommended as particularly suitable, without mentioning or excluding others, that a certain course of action is preferred but not necessarily required, or that (in the negative form) a certain possibility or course of action is undesirable but not prohibited; “may” is used to indicate that a course of action is permissible within the limits of the technical information report; “can” is used as a statement of possibility and capability; “must” is used only for those situations which cannot be otherwise, as in the example “Monday must follow Sunday.”

Suggestions for improving this technical information report are invited. Comments and suggested revisions should be sent to Technical Programs, AAMI, 4301 N. Fairfax Drive, Suite 301, Arlington, VA 22203-1633.

---

NOTE—This introduction does not contain provisions of the AAMI TIR, *Sterilization of health care products—Radiation—Guidance on sterilization of biologics and tissue-based products* (AAMI TIR37:2013), but it does provide important information about the development and intended use of the document.

---

## Introduction

The sources, types and radiation resistances of some microorganisms associated with biologics/tissues might be unique relative to the microorganisms associated with other health care products, such as medical devices. This TIR provides guidance to address issues that are unique to the radiation sterilization of biologics/tissues and to assist in substantiating a sterile label claim for these products. Except where indicated in this TIR, the requirements in the ANSI/AAMI/ISO 11137 series, ANSI/AAMI/ISO 11737-1, and ANSI/AAMI/ISO 11737-2, shall be followed, and AAMI TIR29 is guidance for characterizing the irradiation process and for establishing requisite process controls to ensure the irradiation system remains in a validated state.

The guidance contained in this TIR is only useful if the processing and other steps up to but not including sterilization have already been validated. Examples of steps that should be validated are microbiological testing, cleaning and disinfection, and sterilization of equipment used during processing, packaging and storage.

The key additions/changes in this second edition are as follows:

**1. Change of scope to “biologics and tissue-based products”**

The scope of the document was changed from “human tissue-based products” to “biologics and tissue-based products” because much of the information in the document applies to all types of tissue-based products – not just human – as well as many biological products.

**2. Addition of information for selecting a dose establishment method (Annex A)**

The most significant change to the document was the addition of Annex A, which gives guidance on how to select a dose establishment method that is suitable for a biologic or tissue-based product. These types of products have unique issues, such as availability of samples and tolerance to radiation; therefore, the Annex provides information about certain aspects of the various dose establishment methods and how these aspects might apply to the circumstances associated with a particular product. Some examples of the aspects discussed are the number of products required, availability of batches, bioburden levels versus sterilization dose, and the use of SIPs.

**3. Expansion of information on SIP use in dose establishment**

The concept of SIP was originally developed for traditional medical devices to address issues of size and expected bioburden levels. Most biologics and tissue-based products not only can be tested in their entirety because of their size, but also do not have significant bioburden levels because they are processed in a way that reduces the bioburden. The use of SIPs is not recommended for these products, but guidance is given on what needs to be addressed if an SIP is used for testing.

**4. Addition of a clause on adopting a new product into a family**

Because biologics and tissue-based products are processed in a way that results in very low bioburden, the product family concept is definitely applicable. Although other sterilization methods, such as ethylene oxide, address adopting a new product into a product family, none of the dose establishment documents have guidance about this. A new clause was added that gives guidance on how to assess a new product and show that it can be added to an existing family.

**5. Expansion of information about MPN testing**

Many times biologics and tissue-based products are good candidates for most probable number (MPN) testing, due to their low bioburden and ability to be tested/diluted by weight or volume. Additionally, MPN testing can result in a more sensitive bioburden test with much lower bioburden levels, which allows very low sterilization doses to be established. For these reasons, the information about using MPN testing was expanded upon.

**6. Addition of information on loading patterns for irradiation**

In the new section 5.6.2 (previously 5.5.2), information was added about loading patterns for irradiation. Guidance was included on using a different container for shipment versus irradiation.



# Sterilization of health care products—Radiation— Guidance on sterilization of biologics and tissue- based products

## 1 Scope

### 1.1 Inclusions

This technical information report (TIR) provides guidance for development, validation and routine control associated with the radiation sterilization processing of biologics and tissue-based products.

**NOTE** Although the scope of this document is limited to human biologics and tissue-based products, it provides guidance that might be applicable to other products.

### 1.2 Exclusions

This TIR does not address validation requirements for eliminating and/or inactivating viruses and prions or sterilization of cell/tissue by-products. Guidance on inactivating viruses and prions can be found in ANSI/AAMI/ISO 22442-3:2007. This TIR also does not address handling or safety issues for cell/tissue by-products.

## 2 Normative references

The following normative documents contain provisions which, through reference in this text, constitute provisions of this TIR. For dated references, subsequent amendments to, or revisions of, any of these publications do not apply. However, parties to agreements based on this TIR are encouraged to investigate the possibility of applying the most recent editions of the normative documents indicated below. For undated references, the latest edition of the normative document referred to applies. AAMI maintains a register of currently valid AAMI technical documents.

ANSI/AAMI/ISO 11137-1, *Sterilization of health care products—Radiation—Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices*

ANSI/AAMI/ISO 11137-2, *Sterilization of health care products—Radiation—Part 2: Establishing the sterilization dose*

ANSI/AAMI/ISO 11137-3, *Sterilization of health care products—Radiation—Part 3: Guidance on dosimetric aspects*

ANSI/AAMI/ISO 11737-1, *Sterilization of medical devices—Microbiological methods—Part 1: Determination of a population of microorganisms on products*

ANSI/AAMI/ISO 11737-2, *Sterilization of medical devices—Microbiological methods—Part 2: Tests of sterility performed in the validation of a sterilization process*

This document is intended to be used in conjunction with ANSI/AAMI/ISO 11137-1 and 11137-2.

## 3 Terms and definitions

For the purposes of this document, the following terms and definitions apply.

### 3.1

#### **bacteriostasis/fungistasis (B/F) test (also called the method suitability test or sterility test validation)**

test performed with selected microorganisms to demonstrate the presence or absence of substances that inhibit the multiplication of microorganisms.

**NOTE** For the purposes of this type of testing, the B/F test is considered the validation for the test of sterility and MPN tests.

### **3.2**

#### **batch**

defined quantity of product, intended or purported to be uniform in character and quality, that has been produced during a defined cycle of manufacture.

[ANSI/AAMI/ISO TIR11139:2006, definition 2.1]

NOTE For the purposes of this document, a group of tissue products can be described using terms such as donor, processing batch and sterilization batch.

### **3.3**

#### **bioburden**

population of microorganisms on or in product and/or sterile barrier system.

[ANSI/AAMI/ISO TIR11139:2006, definition 2.2]

### **3.4**

#### **biologics**

a preparation that is synthesized from living organisms, or their products, and used as a diagnostic, preventive or therapeutic agent.

### **3.5**

#### **companion tissue**

tissue from the same donor(s) that is not intended to be used for transplantation.

NOTE For the purposes of this document, companion tissue should be processed in the same manner as tissue that is used for transplantation. Companion tissue is representative of tissue intended for transplantation but is only used for evaluation and/or testing purposes.

### **3.6**

#### **culture conditions**

combination of growth media and manner of incubation used to promote germination, growth and/or multiplication of microorganisms

NOTE The manner of incubation may include the temperature, time and any other conditions specified for incubation.

[ANSI/AAMI/ISO TIR11139:2006, definition 2.10]

### **3.7**

#### **donor identification**

unique code assigned to all tissue and companion tissue that originates from the same donor.

### **3.8**

#### **human tissue-based products**

##### **human cells, tissues, or cellular or tissue-based products (HCT/Ps)**

articles containing or consisting of human cells or tissues that are intended for implantation, transplantation, infusion, or transfer into a human recipient. Examples of HCT/Ps include, but are not limited to, bone, ligament, skin, dura mater, heart valve, cornea, hematopoietic stem/progenitor cells derived from peripheral and cord blood, manipulated autologous chondrocytes, epithelial cells on a synthetic matrix, and semen or other reproductive tissue.

[FDA 21 CFR 1271.3(d)]

### **3.9**

#### **processing**

any activity performed on an HCT/P other than recovery, donor screening, donor testing, storage, labeling, packaging, or distribution, such as testing for microorganisms, preparation, sterilization, steps to inactivate or remove adventitious agents, preservation for storage, and removal from storage.

[FDA 21 CFR 1271.3(ff)]

NOTE For the purposes of this document, cellular-based products are not included.

### **3.10**

#### **processing batch**

product that has been processed in the same manner and at the same time.

NOTE for HCT/Ps, the processing batch is from a single donor source.